

MSG 205 - RED FD14 BIOPACK CONNECT 6 PROCEDURE

6 CONNECT TRANSFER 4 (Duration: ~10 min)

6.1 Notify MCC (Start of step)

Record MET: _____/_____:_____

CMP 6.2 ✓ Air intakes free of debris

6.3 CENTRIF 2 – OFF
CENTRIF 3 – OFF

Incubator 6.4 Transfer 8 Containers to PTCU 4

Cent 3 **CONNECT 1** to PTCU 4
301 0–g

302 **CONNECT 2** to PTCU 4
0–g

Rack **CONNECT 3** to PTCU 4
405 0–g

501 **CONNECT 4** to PTCU 4
0–g

Cent 2 **CONNECT 5** to PTCU 4
201 1–g

203 **CONNECT 6** to PTCU 4
1–g

202 **CONNECT 7** to PTCU 4
1–g

204 **CONNECT 8** to PTCU 4
1–g

Incubator 6.5 Close, latch door

6.6 Notify MCC (End of step)

Record MET:
_____/_____:_____