

MSG 235 - RED FD16 BDS-05 CUE CARD UPDATE

Replace DEORBIT PREPARATION CUE CARD (FRONT) (BDS-05 EOC, CUE CARDS) with the steps below. Cut along dashed lines and tape to BDS-05 Cue Card.

CONNECT PGSC

- Shared PGSC
1. Connect PGSC to ECC, double-click BDS-05 icon in 'Desktop\Shuttle Apps\Spacehab\BDS-05'
 2. Select 'Tasks' / 'Connect to ECC', √ECC detected (MP-19)
 3. √'Fluid Loop Configuration' – 'Perfuse', (otherwise wait 3–10 min)

CHANGE SPIN AND VESSEL WALL RATE

4. Select 'Tasks' / 'ECC Direct Access'
Type "**Hd,1,1**", press [Enter]
√'Event Log' – 'Science Table Override'
 5. Select 'Parameters' / 'Spin Filter Rate'
Type "**3**" for spin filter rate in rpm, select 'OK'
√'Event Log' – 'Spin Filter = 3.000 RPM'
√'Spin Filter Rate' = '3.000' in 'Experiment Parameter' block
 6. Select 'Parameters' / 'Wall Rate'
Type "**0**" for wall rate in rpm, select 'OK'
√'Event Log' – 'Vessel Wall Rate = 0 RPM'
√'Vessel Wall Rate' = '0' in 'Experiment Parameter' block
 7. In BE Viewport Window:
√Spin Filter rotating at approx. 3 RPM
√Vessel not rotating
 8. Notify MCC step 7 complete and data below:
Log 'Perfusion Rate': _____
Log 'Infusion Rate': _____
- AC14
- Shared PGSC
- ### **SHUT DOWN & DISCONNECT PGSC, STOW & NOTIFY MCC**
9. Select 'File' / 'Exit'
 10. Shut down and disconnect PGSC, secure null modem cable, temp stow PGSC
 11. Stow Cue Cards in BDS-05 EOC (CUE CARDS), pg 6–3 through 6–15
 12. Notify MCC (DEORBIT PREPARATION complete)