

## MSG 071 - BLUE FD06 MPFE CARD JAM TROUBLESHOOTING

On MS3:

Move Frozen Saliva from MET 004/14:25 to MET 004/14:15

Move Meal from MET 004/14:40 to MET 004/14:30

Perform the following 5 steps at MET 004/15:30 to MET 004/15:55

(AS07) 1. Perform MP-02 ERROR MESSAGES ON MPFE LECS (JSC-HLS EOC, pg 2-6) up to "Remove and Reinsert Card Tray" in Block 7

(AS07,AC09) 2. Swap Card Tray F4 located in the AMS with Card Tray FC located in the TEHM

(AS07) 3. Finish Block 7 of MP-02 up to "Did error occur within the first 20 minutes?"

4. Perform steps 18-20, 23-30 of SETUP LECS AND INITIALIZE (JSC-HLS EOC, pg 1-5) with AMS Incubator door open. While performing steps 28-29, observe movement of the card tray inside the AMS with Incubator door open.

5. Following step 30, report observations to MCC