

MSG 116A - BLUE FD09 BIOPACK INTERNAL AIR DUCT CLEANING & UPDATES

BIOPACK INTERNAL AIR DUCT FILTER CLEANING

1. cb MAIN POWER - OFF (down)
 2. Cover Air Outlet Filter (bottom side of BPK) completely with Gray Tape
 3. Open freezer/cooler drawer door and pull drawer very slowly out until max out position
 4. Push drawer quickly back in completely
 5. Wait 10 seconds
 6. Repeat steps 3-5 five times
 7. Check the freezer/cooler drawer closed
 8. Remove tape from the air outlet
 9. Check that the air inlet and outlet filters are free of debris
 10. cb MAIN POWER - ON (up) with a mark
-

UPDATE: In BIOPACK EOC, update the following steps:

CONNECT 4: step 4.3, pg 4-5, should read: "Clean air inlet filter with gray tape"

STROMA 4: step 4.2, pg 5-4, should read: "Clean air inlet filter with gray tape"

BACTER 4: step 4.2, pg 6-4, should read: "Clean air inlet filter with gray tape"

BACTER 6: step 6.2, pg 6-7, should read: "Clean air inlet filter with gray tape"

YSTRES 17: step 17.2, pg 7-15, should read: "Clean air inlet filter with gray tape"

YSTRES 26: step 26.2, pg 7-26, should read: "Clean air inlet filter with gray tape"

BIOKIN 4: step 4.2, pg 8-4, should read: "Clean air inlet filter with gray tape"