

## MSG 206A - RED FD14 BDS-05 CUE CARD UPDATE

1 Make the following changes to MEDIA SAMPLE AND PCBA ANALYSIS cue card (BDS-05  
2 Experiment Ops Checklist, pg 6-7):

3  
4 1. Add to unstow from Caddy 1, step 1: Redcaps (2), Mask, BZK/Alcohol wipe

5  
6 2. Replace the bottom half of the back side of the cue card with page 2 of this message.  
7 (Note: the only changes are in the STORE MEDIA SAMPLE section.) You can cut along the  
8 dashed lines and tape it to the BDS-05 Cue Card.

### 9 10 **Rationale:**

11 During BDS Setup the red cap on the outlet port was replaced with a septum for sampling  
12 purposes. The septum needs to be replaced with a red cap to prevent the more brittle  
13 septum from being damaged during reentry and ground handling.

## MSG 206A - RED FD14 BDS-05 CUE CARD UPDATE

### **STORE MEDIA SAMPLE**

- 36. Label Media Sample Bag
- 37. Attach Syringe to media sample bag, store sample, discard Syringe
- 37.1 Don Mask
- 37.2 **Using sterile technique**, replace septum on Outlet Port with Red Cap
- AC14 38. Swab OUTLET PORT and Media Sample Bag Septum with BZK/Alcohol Wipe
- 39. Doff Gloves and Mask, Clean up
- Caddy 1 40. Stow Media Sample Bag in Media Sample Storage Bag

### **LOG PCBA ANALYSIS DATA ON PGSC AND LOG DATA FROM PGSC**

- Shared 41. Select 'Tasks' / 'Log PCBA Data'
- PGSC 42. Enter values recently logged in TABLE 8 (CUE CARDS), pg 6–9 on PGSC  
(Enter "0" for 'Cl', 'BUN' if no values recently logged)  
Verify values entered and select 'OK'
- 43. Log rates in TABLE 8 (CUE CARDS), pg 6–9 using PGSC log rates

### **NOTIFY MCC**

- 44. Notify MCC (data from TABLEs 8 & 9 (CUE CARDS), pg 6–9)

### **CLEAN AIR INTAKE SCREENS**

#### **WARNING**

BE is classed as a Module penetrator and must be crew-tended anytime locker door is not secured closed. Otherwise, it must be pushed back into locker and locker door secured.

- 45. If steps 46–50 were performed from step 8, go to step 51
- AC14 46. ✓Air Intake Screens (2) at bottom and lower left of BE for debris  
Clean using Gray Tape if needed
- AC15 47. Pull out Air Intake Screens (2) on left and right side of ECC
- 48. ✓Screens (2) for debris, clean using Gray Tape if needed
- 49. Reinstall Screens (2) on left and right side of ECC
- 50. If directed to perform steps 46–50 from step 8, go to step 9

### **SHUT DOWN AND DISCONNECT PGSC**

- 51. Select 'File' / 'Exit'
- 52. Shut down and disconnect PGSC, secure null modem cable
- Caddy 1 53. Clean up and stow Trash Caddy