

MSG 011 - RED FD02 MGM ACTIVATION PROCEDURE UPDATES

The change of test cell selection on the MGM UI could occur in three places in the MGM EOC: Step 15 in Activation, Step 14 in AP-02, or Step 15 in AP-03. The procedural change would be the same for all three, as shown below.

EOC Procedure Changes for MGM: Red FD02 and subsequent activities

Crew: MS2

MET: 01/02:45

Change: Activation (MGM, NOM PROC) Step 15, pg. 1-9 to read as follows:

Shared PGSC 15. MGM
Select 'Setup' / 'Select Test Cell'
Select Test Cell
Select test cell letter 'A', 'B', 'C' or 'D'
If test cell A or C, select test cell letter 'A'
If test cell D, select test cell letter 'D'
Select 'OK'
Record Load Cell value 'Load' from PGSC screen in table below

Change AP-02 Step to read:

AP-02 SYSTEM RECOVERY (CRUSH MOTOR POSITION KNOWN), pg. 2-10

14. MGM
Select 'Setup' / 'Select Test Cell'
Select Test Cell
Select test cell letter 'A', 'B', 'C' or 'D'
If test cell A or C, select test cell letter 'A'
If test cell D, select test cell letter 'D'
Select 'OK'
Record Run, Load Cell value 'Load' (from PGSC screen) and MET in table below

Change AP-02 Step to read:

AP-03 ACTIVATION WITH FAILED PRESSURE SENSOR, pg. 2-14

Shared PGSC 15. MGM
Select 'Setup' / 'Select Test Cell'
Select Test Cell
Select test cell letter 'A', 'B', 'C' or 'D'
If test cell A or C, select test cell letter 'A'
If test cell D, select test cell letter 'D'
Select 'OK'
Record Load Cell value 'Load' from PGSC screen in table below